

GLOBACID[®]AF med

ALCOHOL-BASED RAPID DISINFECTANT *ready-to-use solution*

GLOBACID[®]AF med is a ready-to-use disinfectant for rapid disinfection, suitable for use on alcohol resistant surfaces of medical devices. Solution doesn't leave residues on surfaces, thus it is suitable also for the disinfection of surfaces that come in contact with food and food products.

Conforms to requirements in healthcare, food industry, beauty services etc.

PRODUCT PROPERTIES

- fast exposure time
- alcohol-based
- doesn't contain aldehydes, quaternary ammonium compounds and chlorine ions
- doesn't leave residues
- particularly suitable for hard-to-reach places
- doesn't need to be washed off
- doesn't change the visual appearance of the surface

USAGE

Spray the solution or apply with a cloth onto the surfaces, allow to take effect. Surfaces must be completely covered with the solution. Close the cap to prevent the product from evaporation. If the surface tolerance is unknown, try it on an inconspicuous area first. The disinfectant GLOBACID[®]AF med is compatible with most medical devices. Items made of organic glass are not suitable with this product.

EXPOSURE TIME

60 seconds disinfection with mechanical action (wiping)
3 min disinfection with without mechanical action

MICROBIOLOGICAL SPECTRUM

Bactericidal EN 13727	60 seconds
Bactericidal and fungicidal EN 13697	60 seconds
Bactericidal and fungicidal (wipe method) EN 16615	60 seconds
Bactericidal and fungicidal EN 17387	3 minutes
Yeasticidal EN 17387	60 seconds
Fungicidal EN 13624	60 seconds
Tuberculocidal EN 14348	30 seconds
Mycobactericidal EN 14348	60 seconds
Virucidal EN 14476 all enveloped viruses and rotavirus	30 seconds
Virucidal EN 14476 full virucidal	60 seconds
Virucidal EN 16777 full virucidal	60 seconds

CONFORMITY TO STANDARDS

EN 13624, EN 13697, EN 13727, EN 14348, EN 14476, EN 16777, EN 16615, EN 17387, DVV/RKI

STORAGE

See the best before date on the package. Store at temperatures 5-25°C and keep away from direct sunlight to assure the quality of the product.

UTILIZATION

Product is ecologically safe, all active ingredients are biodegradable in flush water. Plastic packaging materials are recyclable. For collection, the plastic packaging is suitable when rinsed, the bottle and cap separated. Spray with its metal parts is also suitable for plastic collection. The packaging box is sorted for cardboard collection. The compatibility of the packaging and the product is quality assured.

SAFETY INFORMATION

See the product safety data sheet for additional information. If in relation to the use of GLOBACID®AF med a serious deterioration of health occurs, it should be reported to the manufacturer and the competent authority of your country. When in doubt, consult with an advisor on regulatory affairs.

TRAINING ADVICE

Thoroughly read the instruction for use and safety data sheet before using this product. For detailed information refer to the technical data sheet.

COMPOSITION

100g of the product contains active ingredients: propan-2-ol 60%, ethanol 15%

Product bears the CE marking according to Medical Device directive 93/42/EEC



PACKAGES: 500ml spray, 1L flip-cap/spray, 5L can



Goodpoint Chemicals OÜ
Urda tee 2/1, Jälgimäe 76404,
Saku vald, Harjumaa, Estonia
www.goodpointchemicals.com
info@goodpointchemicals.com

