

GLOBACID[®]ID

AMINE-BASED CLEANING AND DISINFECTION AGENT FOR INSTRUMENTS

concentrate

Globacid[®]ID is a cost-effective, wide spectrum agent with highly effective properties for removing organic debris; it is an aldehyde, phenol and chlorine free product for the pre- and final disinfection of instruments and accessories made from metal and artificial materials, as well as for the disinfection of all endoscopes.

Also conforms to requirements in food industry, beauty services etc.

PRODUCT PROPERTIES

- wide spectrum of efficacy
 - suitable for use in instrument and ultrasonic baths as well as endoscope washing machines (recommended to add foam blocker)
 - ready-to-use working solution can be preserved in a closed vessel for 4 weeks
-

USAGE

Make the working solution of clean water and concentrate according to the exposure time needed (look at the table). Disassemble the instruments and submerge them into a lid-equipped tub filled with the working solution. Make sure all the instruments are covered with the solution. Once the exposure time has passed, wash or rinse the instruments under running water (at least tap water quality) or with sterile water, as needed and dry with lint-free towel or disposable tissues (instruments can't be left out to dry in the air). It's not recommended to keep the instruments/accessories in the solution longer than the exposure time.

In order to make a working solution out of Globacid[®]ID, use water that is cool or at room temperature.

Note! Do not immerse dental handpieces and turbines in the solution. Process them with a rapid disinfectant (like Vitasept[®]Wipes med) by wiping their outer surfaces.

Note! Avoid contact with aldehyde-based products! Accordingly, surfaces must be cleaned of aldehyde residues. Devices containing butyl rubber parts should be treated with room temperature solution.

WORKING SOLUTION CONCENTRATIONS AND EXPOSURE TIMES

Grossly soiled environment: bactericidal, tuberculocidal, fungicidal(<i>C.albicans</i>), virucidal efficacy against enveloped viruses (incl. HBV, HIV, HCV, coronavirus, influenza virus) and virucidal efficacy against non-enveloped viruses (incl. rota-, noro-, adeno- and poliovirus).	1,25% working solution (5ml/2L)	15min
Low organic load environment: bactericidal, tuberculocidal, fungicidal (<i>C.albicans</i>), virucidal efficacy against	0,25% working solution (25ml/2L)	15min

enveloped viruses (incl. HBV, HIV, HCV, coronavirus, influenza virus) and virucidal efficacy against non-enveloped viruses (adenovirus).	4% working solution (80ml/2L)	5min
Sporicidal efficacy (<i>Clostridium difficile</i>), clean	4% working solution (80ml/2L)	40 min
Fungicidal efficacy (spores of <i>A.brasiliensis</i>) EN 13624, dirty	2% working solution	5 min
Fungicidal efficacy (spores of <i>A.brasiliensis</i>) EN 14562, dirty	2% working solution	30 min

MICROBIOLOGICAL SPECTRUM

Product has bactericidal, tuberculocidal, fungicidal, virucidal efficacy against enveloped viruses (incl. HBV, HIV, HCV, coronavirus, influenza, vaccinia virus) and virucidal efficacy against non-enveloped viruses (incl. rota-, noro-, adeno-, polio- and polyomavirus) and sporidical efficacy (*Clostridium difficile*).

CONFORMITY TO STANDARDS

EN 13624, EN 13727, EN 14348, EN 14476, EN 14561, EN 14562, EN 17126, EN 17111, DVV/RKI

STORAGE

See the best before date on the package. Store at temperatures 5-25°C, keep away from direct sunlight, avoid freezing to assure the quality of the product.

UTILIZATION

The working solution is ecologically safe, all active ingredients are biodegradable in flush water. Plastic packaging materials are recyclable. For collection, the plastic packaging is suitable when rinsed, the bottle and cap separated. The packaging box is sorted for cardboard collection. The compatibility of the packaging and the product is quality assured.

SAFETY INFORMATION

See the product safety data sheet for additional safety information.

COMPOSITION

100g of the product contains active ingredients: didecyldimethylammonium chloride 9g, N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine 14g

pH 11-12

Product bears the CE marking according to Medical Device directive 93/42/EEC  0537

PACKAGES: 1L, 5L



Goodpoint Chemicals OÜ
Urda tee 2/1, Jälgimäe 76404,
Saku vald, Harjumaa, Estonia
www.goodpointchemicals.com
info@goodpointchemicals.com



Use disinfectants safely. Always read the label and product information before use.

09/2024