

GLOBACID®SF med

AMINE-BASED CLEANING AND DISINFECTION AGENT concentrate

GLOBACID®SF med is used for the cleaning and disinfecting of various water-resistant medical devices (patient chairs, stretchers, bed frames etc) and surfaces (doors, windowsills, cupboards, floors, bed frames, kitchen inventory, switches etc). Suitable for cleaning and disinfecting operating room inventory.

Also conforms to requirements in food industry, beauty services etc.

PRODUCT PROPERTIES

- cost-efficient
- wide spectrum of efficacy
- aldehyde-free
- great cleaning properties
- material friendly (metal, plastic, plexi)
- neutral scent

USAGE

Make the working solution of clean water and concentrate according to the exposure time needed (look at the table). In order to make a working solution and for better cleaning properties, use water that is cool or at room temperature. Once the exposure time has passed, all surfaces that come in contact with food must be rinsed.

NB! Avoid contact with aldehyde-based products. Accordingly, surfaces previously treated with aldehyde-based products must be cleaned from aldehyde residues.

EXPOSURE TIME

- General disinfection, preliminary cleaned environment 0,25% = 15min
- Infection cases, dirty environment 1,25% = 15min

MICROBIOLOGICAL SPECTRUM

| | | | |
|---|-------|----------|--------|
| General disinfection, preliminary cleaned environment incl. tuberculocidal* | 0,25% | 2,5ml/l | 15 min |
| Infection cases, dirty environment, bactericidal, fungicidal (<i>Candida albicans</i>), tuberculocidal, full virucidal** | 1,25% | 12,5ml/l | 15 min |
| Fungicidal (<i>C. albicans</i>, <i>A. brasiliensis</i>)** | 2,5% | 25ml/l | 5 min |
| Bactericidal, fungicidal (<i>Candida albicans</i>), all enveloped viruses (HBV, HCV, HIV)*, adeno virus*, yeasticidal (<i>Candida albicans</i>)* | 0,25% | 2,5ml/l | 5 min |
| Sporicidal efficacy (<i>Clostridium difficile</i>)* | 2,5% | 25ml/l | 15 min |

*clean conditions

**dirty conditions

CONFORMITY TO STANDARDS

EN 1276, EN 13624, EN 13727, EN 14348, EN 14476, EN 14561, EN 14562, EN 16615, EN 17387, EN 17126, DVG, DVV/RKI, DGHM/ÖGHMP

STORAGE

See the best before date on the package. Store at temperatures 5-25°C, keep away from direct sunlight, avoid freezing to assure the quality of the product.

UTILIZATION

The working solution is ecologically safe, all active ingredients are biodegradable in flush water. Plastic packaging materials are recyclable. For collection, the plastic packaging is suitable when rinsed, the bottle and cap separated. The packaging box is sorted for cardboard collection. The compatibility of the packaging and the product is quality assured.

SAFETY INFORMATION

See the product safety data sheet for additional safety information.

COMPOSITION

100g of the product contains active ingredients: didecyldimethylammonium chloride 6g, N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine 5,5g

pH 10 - 11

Product bears the CE marking according to Medical Device directive 93/42/EEC
Biocide reg.no: 2251/21



PACKAGES: 1L dosage bottle, 5L canister, 6L canister, 10L canister



Goodpoint Chemicals OÜ
Urda tee 2/1, Jälgimäe 76404,
Saku vald, Harjumaa, Estonia
www.goodpointchemicals.com
info@goodpointchemicals.com

