

GLOBATABS®

CHLORINE TABLETS

GLOBATABS® is used for disinfecting surfaces, tools and equipment in health care, social welfare, nursing and educational institutions, saunas, pools and the food industry; for rendering harmless surfaces contaminated with blood and other bodily fluids.

PRODUCT PROPERTIES

- accurate measured dose of stabile active chlorine is released when dissolved in water
 - broad biocidal spectrum against bacteria, fungi, viruses and spores
 - research shows that sodium dichloroisocyanurate dihydrate solutions are active and effective up to 10 weeks when stored in closed, clean container and kept away from direct sunlight!
 - several times more effective compared to regular hypochlorite containing products
 - very good cleaning properties thanks to the sodium bicarbonate contained in the product
-

USAGE

Surface disinfection and unpleasant scent removal

Add 1 tablet to 7 litres of water (200 ppm). Apply on surfaces with a mop, sponge or cloth. Allow 10 minutes contact time and rinse surfaces with fresh water.

* Follow the local regulation to use correct dosage of active ingredient in ppm

Soiled mops disinfection

Add 1 tablet to 0,75 litres of water (2000 ppm). Allow 15 minutes contact time and rinse mops with fresh water.

* Follow the local regulation to use correct dosage of active ingredient in ppm

Sanitation of egg shells

Add 1 tablet to 7 litres of water (200 ppm). Dip eggs in the solution and remove from solution. Allow 10 minutes contact time and rinse with drinking water.

* Follow the local regulation to use correct dosage of active ingredient in ppm

Sanitation of salad/vegetables

Add 1 tablet to 7 litres of water (200 ppm). Ensure the tablet is fully dissolved before adding salad/vegetables. Allow 10 minutes contact time and rinse thoroughly with drinking water.

* Follow the local regulation to use correct dosage of active ingredient in ppm

Disposal of used syringes, needles etc

Add 1 tablet to 0,75 litres of water (2000 ppm). Pour the solution in the needle and syringe disposal unit.

* Follow the local regulation to use correct dosage of active ingredient in ppm

Disposal of medical waste (tissues, blood etc)

Add 2 tablets to 075 litres of water (4000 ppm).

* Follow the local regulation to use correct dosage of active ingredient in ppm

Active/free chlorine in solution		Tablets (pcs)	Water (L)
ppm	%		
200	0,02	1	7
1000	0,1	1	1,5
2000	0,2	1	0,75
4000	0,4	2	0,75

Swimming pool chlorination

Active chlorine is the part of chlorine present in water that hasn't reacted with other organic materials, it's added in the water and has a great disinfecting property, to keep the water clean and clear. Active chlorine doesn't cause skin nor eye irritation, the amount needed depends on water temperature, sunshine, how many swimmers are in the water and how contaminated the pool is. For ensuring pool water safety for swimmers, the water should be disinfected regularly so the active chlorine content is 0,5...2 ppm (mg/l). Water pH, active and inactive chlorine amount, cloudiness and temperature should be measured every day, several times a day. Follow the local regulation to use correct dosage of active ingredient in ppm.

Check table for dosing.

Active/free chlorine in solution		Tablets (pcs)	Water (L)
ppm	%		
0,5	0,00005	1	3000
1	0,0001	1	1500
1,5	0,00015	1	1000
2	0,0002	2	1500

EXPOSURE TIME

Bactericidal, fungicidal, virucidal in medical dirty conditions	0,2%	2000ppm	5 min
Mycobactericidal in medical dirty conditions	0,2%	2000ppm	15 min
Sporicidal efficacy (<i>Clostridioides difficile</i>) in medical dirty conditions	0,4%	4000ppm	15 min
Sporicidal efficacy (<i>B.subtilis</i> , <i>B. cereus</i>) in medical dirty conditions	0,4%	4000ppm	30 min

CONFORMITY TO STANDARDS

EN 13624, EN 13727, EN 16615 (60 sec), EN 17387, EN 14348, EN 14476, EN 17126

STORAGE

See the best before date on the package. Store at temperatures 5-25°C, in a well-ventilated room, away from heat and ignition sources and direct sunlight.

UTILIZATION

The working solution is ecologically safe, all active ingredients are biodegradable in flush water. Plastic packaging materials are recyclable. For collection, the plastic packaging is suitable when rinsed, the bottle and cap separated. The packaging box is sorted for cardboard collection. The compatibility of the packaging and the product is quality assured.

SAFETY INFORMATION

See the product safety data sheet for additional safety information.

COMPOSITION


1 tablet contains: Sodium dichloroisocyanurate dihydrate (CAS: 51580-86-0) 2,7g (1,5g active chlorine), sodium bicarbonate, adipic acid.

pH ca 6

Use disinfectants safely. Always read the label and product information before use.

04/2024

PACKAGES



300 tablets in a jar



Goodpoint Chemicals OÜ
Urda tee 2/1, Jälgimäe 76404,
Saku vald, Harjumaa, Estonia
www.goodpointchemicals.com
info@goodpointchemicals.com