


Directions for use

REF: 603.003.0001

Sterintech™ 3 - line Label gun

Content:

- 1) Installing the ink roller
- 2) Removing the ink roller
- 3) Loading the labels
- 4) Setting print characters
- 5) Clearing label jams, cleaning and lubrication
- 6) Lay-out of the 3 lines digits



1) Installing the ink roller

Open the label gun cover by pressing the knob at the top front of the label gun.

Open the package in which the ink roller is packed and take out the ink roller

Hold the new ink roller on each end in between your finger tips

Snap the ink roller into place.

2) Removing the ink roller

Open the label gun cover by pressing the knob at the top front of the label gun

Hold the ink roller on each end and push the release bar to release the ink roller


Discard the used, old ink roller. As it is non-toxic the old ink roller can be treated as normal waste.

3) Loading labels

Open the label gun cover by pressing the knob at the top front of the label gun. Swing the cover completely backwards until you hear a click. Don't force the cover past the click position.

Remove the closing label from the label roll.

Peel and discard the first 10 labels from the backing paper.



cont'd : 3) Loading labels.

Place the label roll between the supply hubs. Do not thread label strip through any rollers.

Drape the supply across the label gun and center it over the peel bar. The supply should hang out about 15 cm past the peel bar in order to comfortably feed in the labels.

Close the cover

Feed the label strip under the large rollers and into the slot above the arrows.

Hold the strip in place and squeeze the trigger until the label strip feeds through the rear exit.


4) Setting print characters

Pull the knob out until the character pointers are on the desired position.

Turn the knob until the desired character with white background appears between the pointers

Repeat these two steps for each position

Push in the knob when you are finished.



5) Clearing label jams, cleaning and lubrication


Open the cover and remove the label roll

Remove any loose labels from the inside of the label gun

Pull up the deflector frame to expose the backing paper path

Remove any loose labels from under the deflector frame

Snap back the deflector frame back into place.



Use a commercial adhesive solvent to clean adhesive build-up

Reload and close the cover.

In case needed you may apply 1 drop of 100% silicone lubricant on the places marked by the arrows in the picture at the left.

6) Lay-out of the 3 lines digits

TOP LINE												MIDDLE LINE												BOTTOM LINE														
Pos.	1	2	3	4	5	6	7	8	9	10	11	12	Pos.	1	2	3	4	5	6	7	8	9	10	11	12	Pos.	1	2	3	4	5	6	7	8	9	10	11	12
REF: 603.003.0001																																						
Steritech 3 - line Label gun (12+12+12 digits)																																						
standard lay-out																																						
1													1													1												
2	A	A	A	A	A	A	A	A	2	A	A	A	A	A	A	A	A	2
3	1	1	1	1	B	B	B	B	B	B	B	B	3	1	1	1	1	B	B	B	B	B	B	B	B	3	1	1	1	1	1	1	1	1	1	1	1	1
4	2	2	2	2	C	C	C	C	C	C	C	C	4	2	2	2	2	C	C	C	C	C	C	C	C	4	2	2	2	2	2	2	2	2	2	2	2	2
5	3	3	3	3	D	D	D	D	D	D	D	D	5	3	3	3	3	D	D	D	D	D	D	D	D	5	3	3	3	3	3	3	3	3	3	3	3	3
6	4	4	4	4	E	E	E	E	E	E	E	E	6	4	4	4	4	E	E	E	E	E	E	E	E	6	4	4	4	4	4	4	4	4	4	4	4	4
7	5	5	5	5	F	F	F	F	F	F	F	F	7	5	5	5	5	F	F	F	F	F	F	F	F	7	5	5	5	5	5	5	5	5	5	5	5	5
8	6	6	6	6	G	G	G	G	G	G	G	G	8	6	6	6	6	G	G	G	G	G	G	G	G	8	6	6	6	6	6	6	6	6	6	6	6	6
9	7	7	7	7	H	H	H	H	H	H	H	H	9	7	7	7	7	H	H	H	H	H	H	H	H	9	7	7	7	7	7	7	7	7	7	7	7	7
10	8	8	8	8	J	J	J	J	J	J	J	J	10	8	8	8	8	J	J	J	J	J	J	J	J	10	8	8	8	8	8	8	8	8	8	8	8	8
11	9	9	9	9	K	K	K	K	K	K	K	K	11	9	9	9	9	K	K	K	K	K	K	K	K	11	9	9	9	9	9	9	9	9	9	9	9	9
12	0	0	0	0	L	L	L	L	L	L	L	L	12	0	0	0	0	L	L	L	L	L	L	L	L	12	0	0	0	0	0	0	0	0	0	0	0	0
13	M	M	M	M	M	M	M	M	M	M	M	M	13	M	M	M	M	M	M	M	M	M	M	M	M	13	FOR											
14	C	C	C	C	N	N	N	N	N	N	N	N	14	C	C	C	C	N	N	N	N	N	N	N	N	14	\$	\$	\$	\$	\$	\$	\$	\$	\$	\$	\$	\$
15	G	G	G	G	P	P	P	P	P	P	P	P	15	G	G	G	G	P	P	P	P	P	P	P	P	15	/	/	/	/	/	/	/	/	/	/	/	/
16	A	A	A	A	Q	Q	Q	Q	Q	Q	Q	Q	16	A	A	A	A	Q	Q	Q	Q	Q	Q	Q	Q	16	-	-	-	-	-	-	-	-	-	-	-	-
17	pm	pm	pm	pm	R	R	R	R	R	R	R	R	17	pm	pm	pm	pm	R	R	R	R	R	R	R	R	17												
18	%	%	%	%	S	S	S	S	S	S	S	S	18	%	%	%	%	S	S	S	S	S	S	S	S	18												
19													19													19												
20					T	T	T	T	T	T	T	T	20					T	T	T	T	T	T	T	T	20												
21					U	U	U	U	U	U	U	U	21					U	U	U	U	U	U	U	U	21												
22					V	V	V	V	V	V	V	V	22					V	V	V	V	V	V	V	V	22												
23					W	W	W	W	W	W	W	W	23					W	W	W	W	W	W	W	W	23												
24					X	X	X	X	X	X	X	X	24					X	X	X	X	X	X	X	X	24												
25					Y	Y	Y	Y	Y	Y	Y	Y	25					Y	Y	Y	Y	Y	Y	Y	Y	25												
26					Z	Z	Z	Z	Z	Z	Z	Z	26					Z	Z	Z	Z	Z	Z	Z	Z	26												
27					1	1	1	1	1	1	1	1	27					1	1	1	1	1	1	1	1	27												
28					2	2	2	2	2	2	2	2	28					2	2	2	2	2	2	2	2	28												
29					3	3	3	3	3	3	3	3	29					3	3	3	3	3	3	3	3	29												
30					4	4	4	4	4	4	4	4	30					4	4	4	4	4	4	4	4	30												
31					5	5	5	5	5	5	5	5	31					5	5	5	5	5	5	5	5	31												
32					6	6	6	6	6	6	6	6	32					6	6	6	6	6	6	6	6	32												
33					7	7	7	7	7	7	7	7	33					7	7	7	7	7	7	7	7	33												
34					8	8	8	8	8	8	8	8	34					8	8	8	8	8	8	8	8	34												
35					9	9	9	9	9	9	9	9	35					9	9	9	9	9	9	9	9	35												
36					0	0	0	0	0	0	0	0	36					0	0	0	0	0	0	0	0	36												