

CYCLONE RTU med

DISINFECTANT WITH CLEANING PROPERTIES

ready-to-use solution

Cyclone RTU med is an effective product for cleaning medical devices (bedpans, commode seats etc) and for removing unpleasant odour. The product is suitable for cleaning water-resistant surfaces.

Conforms to requirements in healthcare, food industry, beauty services etc.

PRODUCT PROPERTIES

- good cleaning properties
 - removes unpleasant odour, even from urine
 - rapid effect (1 minute)
-

USAGE

Spray the solution or apply with a cloth onto the surfaces, allow to take effect for 1 minute and rinse with clean water.

EXPOSURE TIME

60 seconds

MICROBIOLOGICAL SPECTRUM

MICROORGANISM	SOILING CONDITIONS	CONCENTRATION	CONTACT TIME
Bactericidal	dirty conditions	undiluted	60 seconds
Yeasticidal			
<i>Candida albicans</i>	dirty conditions	undiluted	60 seconds

CONFORMITY TO STANDARDS

EN 13727, EN 1275, EN 13624

STORAGE

See the best before date on the package. Store at room temperature, keep away from direct sunlight and avoid freezing to assure the quality of the product.

UTILIZATION

Product is ecologically safe, all active ingredients are biodegradable in flush water. Plastic packaging materials are recyclable. For collection, the plastic packaging is suitable when rinsed, the bottle and cap separated. Spray

with its metal parts is also suitable for plastic collection. The packaging box is sorted for cardboard collection. The compatibility of the packaging and the product is quality assured.

SAFETY INFORMATION

See the product safety data sheet for additional safety information.

COMPOSITION

100g of the product contains active ingredients: benzyl-C12-18-alkyldimethylchloride 0,12g, C12-14-alkyl[(ethylphenyl)methyl]dimethylchloride 0,12g

pH 3

Product bears the CE marking according to Medical Device directive 93/42/EEC




PACKAGES



1L spray bottle



1L flip-cap bottle



5L canister



Goodpoint Chemicals OÜ
Urda tee 2/1, Jälgimäe 76404,
Saku vald, Harjumaa, Estonia
www.goodpointchemicals.com
info@goodpointchemicals.com

