

GLOBACID[®]SF med

AMINE-BASED CLEANING AND DISINFECTING WIPES

Used for the cleaning and disinfecting of various water-resistant medical devices (patient chairs, stretchers, bed frames etc) and surfaces (doors, windowsills, cupboards, floors, bed frames, kitchen inventory, switches etc). Suitable for cleaning and disinfecting operating room inventory.

Also conforms to requirements in food industry, beauty services etc.

PRODUCT PROPERTIES

- cost-efficient
- wide spectrum of efficacy
- aldehyde-free
- great cleaning properties
- material friendly (metal, plastic, plexi)
- neutral scent

USAGE

Pull the wipes out until perforated spot according to the amount necessary and if needed support with other hand in order to tear away. All surfaces should be wiped, make sure that the surface is covered with disinfectant. Throw the used wipe away. In order to avoid the drying of wipes, close the flip.

NB! Avoid contact with aldehyde-based products. Accordingly, surfaces previously treated with aldehyde-based products must be cleaned from aldehyde residues.

EXPOSURE TIME

60 seconds bactericidal, yeasticidal (EN 16615) dirty conditions

5 minutes fungicidal (*C.albicans*, *A. brasiliensis*) dirty conditions

15 minutes mycobactericidal (*M. terrae*, *M.avium*) dirty conditions

15 minutes full virucidal dirty conditions

15 minutes sporicidal (*Clostridium difficile*) efficacy, clean conditions

MICROBIOLOGICAL SPECTRUM

Product has bactericidal, mycobactericidal, fungicidal, virucidal efficacy against enveloped viruses (incl. HBV, HIV, HCV, coronavirus, influenza, vaccinia virus) and virucidal efficacy against non-enveloped viruses (incl. rota-, noro-, adeno-, polio- and polyomavirus) and sporicidal efficacy (*Clostridium difficile*).

CONFORMITY TO STANDARDS

EN 1276, EN 13624, EN 13697, EN 13727, EN 14348, EN 14476, EN 14561, EN 14562, EN 16615, EN 17126

STORAGE

See the best before date on the package. Store at temperatures 5-25°C, keep away from direct sunlight, avoid freezing to assure the quality of the product.

UTILIZATION

Plastic packaging materials are recyclable. For recycling, rinse package with water. The packaging box is sorted for cardboard collection. The compatibility of the packaging and the product is quality assured.

SAFETY INFORMATION

See the product safety data sheet for additional safety information.

COMPOSITION

100g of the product contains active ingredients: didecyldimethylammonium chloride 0,3g, N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine 0,28g

pH 8

Product bears the CE marking according to Medical Device directive 93/42/EEC  0537

PACKAGES: 24pcs flow bag, 120pcs flow bag/jar, 150pcs jar, 200pcs flow bag/jar, 600pcs jar



Goodpoint Chemicals OÜ
Urda tee 2/1, Jälgimäe 76404,
Saku vald, Harjumaa, Estonia
www.goodpointchemicals.com
info@goodpointchemicals.com

