

GLOBACID®STERIL

CHEMICAL STERILANT FOR INSTRUMENTS *ready-to-use solution*

Globacid[®]Steril is an aldehyde-based, phenol- and chlorine-free agent for high-level disinfection and chemical sterilization of instruments and accessories made of metal and all synthetic materials. Suitable for sterilizing instruments and accessories used in beauty services.

PRODUCT PROPERTIES

- wide spectrum of efficacy
- user-friendly
- depending on the workload, the working solution in use can be used for 30 days

USAGE

Pre-washed, disinfected (e.g. with Globaclean[®]) and disassembled instruments should be placed into a lidequipped tub filled with the working solution. Make sure all the instruments are covered with the solution. Once the exposure time has passed, wash or rinse the instruments under running water (at least tap water quality) or with sterile water, as needed and dry with lint-free towel or disposable tissues (instruments can't be left out to dry in the air). It's not recommended to keep the instruments/accessories in the solution longer than the exposure time.

Note! Avoid contact with amine-based products! Instruments <u>should not be</u> disinfected with Globacid[®]ID before placing them in Globacid[®]Steril solution.

EXPOSURE TIME

Chemical sterilisation - 30 minutes

MICROBIOLOGICAL SPECTRUM

The product has a full sporicidal effect (*C.difficile, B.subtilis, B.cereus*), which means that all microorganisms (including fungi, bacteria, viruses) are destroyed.

CONFORMITY TO STANDARDS

EN 17126

STORAGE

See the best before date on the package. Store at temperatures 5-25°C, keep away from direct sunlight, avoid freezing to assure the quality of the product.

UTILIZATION

Product is ecologically safe, all active ingredients are biodegradable in flush water. Plastic packaging materials are recyclable. For collection, the plastic packaging is suitable when rinsed, the bottle and cap separated. The packaging box is sorted for cardboard collection. The compatibility of the packaging and the product is quality assured.

SAFETY INFORMATION

See the product safety data sheet for additional safety information.

COMPOSITION

100g of the product contains active ingredients: glutaral 2,6g, didecyldimethylammonium chloride 0,85g

pH 5-7

PACKAGES





Goodpoint Chemicals OÜ Urda tee 2/1, Jälgimäe 76404, Saku vald, Harjumaa, Estonia www.goodpointchemicals.com info@goodpointchemicals.com

