

GLOBACID®ID ready-to-use

AMINE-BASED CLEANING AND DISINFECTION AGENT FOR INSTRUMENTS

ready-to-use solution

Globacid®ID is a cost—effective, wide spectrum agent with highly effective properties for removing organic debris; it is an aldehyde, phenol and chlorine free product for the pre— and final disinfection of instruments and accessories made from metal and artificial materials, as well as for the disinfection of all endoscopes.

Also conforms to requirements in food industry, beauty services etc.

PRODUCT PROPERTIES

- wide spectrum of efficacy
- suitable for use in instrument and ultrasonic baths as well as endoscope washing machines (recommended to add foam blocker)

USAGE

Spray the product onto disassembled instruments or submerge them into a lid-equipped tub filled with the working solution. Make sure all the instruments are covered with the solution. Once the exposure time has passed, wash or rinse the instruments under running water (at least tap water quality) or with sterile water, as needed and dry with lint-free towel or disposable tissues (instruments can't be left out to dry in the air). It's not recommended to keep the instruments/accessories in the solution longer than the exposure time.

Note! Do not immerse dental handpieces and turbines in the solution. Process them with a rapid disinfectant (like Vitasept Wipes med) by wiping their outer surfaces.

Note! Avoid contact with aldehyde-based products! Accordingly, surfaces must be cleaned of aldehyde residues. Accordingly, surfaces must be cleaned of aldehyde residues. Devices containing butyl rubber parts should be treated with room temperature solution.

EXPOSURE TIME

5 minutes (Clostridium difficile 40 min)

MICROBIOLOGICAL SPECTRUM

Product has bactericidal, tuberculocidal, fungicidal, virucidal efficacy against enveloped viruses (incl. HBV, HIV, HCV, coronavirus, influenza, vaccinia virus) and virucidal efficacy against non-enveloped viruses (incl. rota-, noro, adeno-, polio- and polyomavirus) and sporidical efficacy (*Clostridium difficile*).

CONFORMITY TO STANDARDS

EN 13624, EN 13727, EN 14348, EN 14476, EN 14561, EN 14562, EN 17126, EN 17111, DVV/RKI

STORAGE

See the best before date on the package. Store at temperatures 5-25°C, keep away from direct sunlight, avoid freezing to assure the quality of the product.

UTILIZATION

Product is ecologically safe, all active ingredients are biodegradable in flush water. Plastic packaging materials are recyclable. For collection, the plastic packaging is suitable when rinsed, the bottle and cap separated. Spray with its metal parts is also suitable for plastic collection. The packaging box is sorted for cardboard collection. The compatibility of the packaging and the product is quality assured.

SAFETY INFORMATION

See the product safety data sheet for additional information. If in relation to the use of Globacid®ID a serious deterioration of health occurs, it should be reported to the manufacturer and the competent authority of your country. When in doubt, consult with an advisor on regulatory affairs.

TRAINING ADVICE

Thoroughly read the instruction for use and safety data sheet before using this product. For detailed information refer to the technical data sheet.

COMPOSITION

100g of the product contains active ingredients: N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine 0,56g, didecyldimethylammonium chloride 0,36g

pH 10-11



PACKAGES: 1L spray, 5L



Goodpoint Chemicals OÜ Urda tee 2/1, Jälgimäe 76404, Saku vald, Harjumaa, Estonia www.goodpointchemicals.com info@goodpointchemicals.com

