

GLOBADENT[®]RTU

ALDEHYDE-FREE AGENT FOR DISINFECTING DENTAL IMPRESSIONS AND DENTURES

ready-to-use solution

Globadent[®]RTU is an amine-based ready-to-use solution for washing and disinfecting dental impressions, dentures, dental bridges, dental crowns, and dental spoons. Used for cleaning and disinfecting silicone, alginate, polyether rubber, polysulfide and hydrocolloid impressions, prosthetic works, etc

PRODUCT PROPERTIES

- wide spectrum of efficacy
 - quick exposure time (5 min)
 - plaster cast model accuracy thanks to exceptional cleaning power (blood, saliva, etc)
-

USAGE

After taking a dental impression or trying the dentures, immerse them for 5 minutes in a container closed with a lid and filled with the working solution. Impressions grossly soiled with blood should be rinsed before immersing in the solution. Make sure that all surfaces that are disinfected are covered with the solution. After the exposure period has elapsed wash the impressions/dentures with water (at least drinking water quality). The solution may be used for up to 7 working days, depending on the workload.

Note! It is not recommended to keep the dental impressions/dentures etc in the solution beyond the exposure time.

Note! Avoid contact with glutaraldehyde!

EXPOSURE TIME

5 minutes

MICROBIOLOGICAL SPECTRUM

Product has bactericidal, tuberculocidal, fungicidal and virucidal efficacy (incl. HBV, HIV, HCV, coronavirus, influenza, vaccinia virus).

CONFORMITY TO STANDARDS

EN 1276, EN 13727, EN 14348, EN 14476, EN 14562, DVV/RKI, DGHM/ÖGHMP

STORAGE

See the best before date on the package. Store at temperatures 5-25°C, cap tightly closed, keep away from direct sunlight, avoid freezing to assure the quality of the product.

UTILIZATION

Product is ecologically safe, all active ingredients are biodegradable in flush water. Plastic packaging materials are recyclable. For collection, the plastic packaging is suitable when rinsed, the bottle and cap separated. The packaging box is sorted for cardboard collection. The compatibility of the packaging and the product is quality assured.

SAFETY INFORMATION

See the product safety data sheet for additional safety information.

COMPOSITION

100g of the product contains active ingredients: N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine 0,56g, didecyldimethylammonium chloride 0,36g

pH 10-11

Product bears the CE marking according to Medical Device directive 93/42/EEC



PACKAGES: 5L can



Goodpoint Chemicals OÜ
Urda tee 2/1, Jälgimäe 76404,
Saku vald, Harjumaa, Estonia
www.goodpointchemicals.com
info@goodpointchemicals.com

