

GLOBAQUAT[®] RTU med

ALCOHOL-FREE CLEANING AND RAPID DISINFECTION AGENT

ready-to-use solution

Used for cleaning and disinfecting acrylic surfaces, different types of plastic, fitness club inventory, cables, toys, mattresses, medical devices (ECG, ultrasound probes, infant incubators, dental patient chairs incl. their cables, hygienic and anti-decubitus mattresses, etc.).

Conforms to requirements in healthcare, food industry, beauty services etc.

PRODUCT PROPERTIES

- alcohol- and aldehyde-free
- odorant- and colorant-free
- quick exposure time (1 min)
- suitable for various kind of materials
- excellent cleaning performance
- convenient spray bottle
- not tuberculocidal

USAGE

Spray the solution or apply with a cloth onto the surfaces, allow to take effect for 1 minute. Surfaces have to be completely covered with the solution.

EXPOSURE TIME

60 seconds

MICROBIOLOGICAL SPECTRUM

MICROORGANISM	SOILING CONDITIONS	CONCENTRATION	CONTACT TIME
Bactericidal	clean conditions	undiluted	60 seconds
Fungicidal			
<i>Candida albicans</i>	clean conditions	undiluted	60 seconds
Virucidal			
HBV/HCV, HIV	clean conditions	undiluted	60 seconds
Herpes simplex	clean conditions	undiluted	60 seconds
Influenza viruses	clean conditions	undiluted	60 seconds
Vaccinia virus	clean conditions	undiluted	60 seconds
Coronaviruses	clean conditions	undiluted	60 seconds
Adenovirus	clean conditions	undiluted	60 seconds
Norovirus	clean conditions	undiluted	30 minutes
Rotavirus	clean conditions	undiluted	60 seconds

CONFORMITY TO STANDARDS

EN 1040, EN 1275, EN 1276, EN 13624, EN 13697, EN 13727, EN 14476, EN 14561, EN 14562, EN 16615, DVG, DVV/RKI, DGHM, „Effect on Acrylic Glass“ Röhm GmbH

STORAGE

See the best before date on the package. Store at room temperature, keep away from direct sunlight and avoid freezing to assure the quality of the product.

UTILIZATION

Product is ecologically safe, all active ingredients are biodegradable in flush water. Plastic packaging materials are recyclable. For collection, the plastic packaging is suitable when rinsed, the bottle and cap separated. Spray with its metal parts are also suitable for plastic collection. The packaging box is sorted for cardboard collection. The compatibility of the packaging and the product is quality assured.

SAFETY INFORMATION

See the product safety data sheet for additional safety information.

COMPOSITION


100g of the product contains active ingredients: didecyldimethylammonium chloride 0,25g; benzyl-C12-16-alkyldimethyl, chlorides 0,25g; C12-14-alkyl[(ethylphenyl)methyl]dimethyl, chlorides 0,25g

pH 7


Product bears the CE marking according to Medical Device directive 93/42/EEC




PACKAGES




250ml foam bottle



1L spray bottle



1L flip-cap bottle



5L canister



Goodpoint Chemicals OÜ
Urda tee 2/1, Jälgimäe 76404,
Saku vald, Harjumaa, Estonia
www.goodpointchemicals.com
info@goodpointchemicals.com



Use disinfectants safely. Always read the label and product information before use.

12/2022