GLOBAQUAT® WIPES med

ALCOHOL-FREE CLEANING AND DISINFECTING WIPES

Medical device · Class I

INTENDED USE AS A MEDICAL DEVICE

GLOBAQUAT®WIPES med is intended for cleaning the external surfaces of non-invasive medical devices in healthcare facilities, such as acrylic and plastic surfaces, cables, external surfaces of ECG and ultrasound probes, incubators, dental patient chairs (including their cables), hygienic and anti-decubitus mattresses and similar devices.

INTENDED USE AS A BIOCIDE

For cleaning and disinfection of acrylic surfaces, various plastic surfaces, fitness equipment, cables, toys, mattresses and other moisture-resistant surfaces in healthcare, the food industry, beauty services and other sectors.

PRODUCT PROPERTIES

- alcohol- and aldehyde-free
- odorant- and colorant-free
- cleans and disinfects in one step (disinfection as biocidal use)
- bactericidal, fungicidal and virucidal (biocidal use)
- rapid disinfection: contact time 1 min (biocidal use)
- soft fabric-like material, contains no cellulose and does not leave lint on surfaces
- suitable for various kinds of materials
- excellent cleaning performance
- not tuberculocidal

DIRECTIONS FOR USE

Pull the wipe from the pack to the perforation according to the required amount and, if necessary, support with the other hand to facilitate tearing.

Cleaning of medical devices: thoroughly wipe the external surfaces of the device with a GLOBAQUAT®WIPES med wipe and allow to air-dry.

Disinfection of surfaces (biocidal use): thoroughly wipe the surface and ensure that the area to be treated is evenly covered with the solution. Allow to act for 1 minute and let air-dry.

Dispose of the wipe after use. To prevent the wipes from drying out, always close the lid after use. Not suitable for cleaning hands. Do not use on human or animal skin.

SPECTRUM OF ACTIVITY (BIOCIDAL USE)

Bactericidal, fungicidal (*Candida albicans*) and virucidal (HBV, HCV, HIV, herpes simplex virus, influenza viruses, vaccinia virus, coronaviruses, adenovirus, rotavirus). Undiluted, under clean conditions, contact time 60 seconds.

CONFORMITY TO STANDARDS (BIOCIDAL EFFICACY)

EN 13624, EN 13727, EN 14476, EN 14561, EN 14562, EN 16615

STORAGE

See the best before date on the package. Store at +5 °C to +25 °C, keep away from direct sunlight and avoid freezing. Close the lid after use.

WARNINGS AND PRECAUTIONS

Avoid contact with eyes and mucous membranes. Do not use on human or animal skin. H412 Harmful to aquatic life with long lasting effects. P102 Keep out of reach of children. P273 Avoid release to the environment. P501 Dispose of contents and/or container in accordance with local regulations. See the safety data sheet for additional safety information.

REPORTING OF SERIOUS INCIDENTS

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and to the competent authority of the country in which the user and/or patient is established. If medical advice is needed, have the product container or label at hand.

UTILIZATION

Plastic packaging materials are recyclable. For recycling, rinse the package with water. The packaging box is sorted for cardboard collection. The compatibility of the packaging and the product is quality assured.

COMPOSITION

100 g of product contains active substances: didecyldimethylammonium chloride 0.25 g; C12–16-alkyldimethylbenzylammonium chloride 0.25 g; (C12–14)-alkyldimethyl(ethylbenzyl)ammonium chloride 0.25 g. pH ~7.

REGULATORY INFORMATION

GLOBAQUAT® WIPES med is a Class I medical device according to Regulation (EU) 2017/745, Annex VIII, Rule 1, for its intended use as a medical device (cleaning). The disinfecting use described above is a biocidal use (biocide reg. no. (EE) 2253/21).

CE MD

PACKAGES

120 wipes (flow bag), 150 wipes in a can, 600 wipes in a bucket



Goodpoint Chemicals OÜ
Urda tee 2/1, Jälgimäe 76404,
Saku vald, Harjumaa, Estonia
www.goodpointchemicals.com
info@goodpointchemicals.com